


**EN** Titanium Forceps

**DE** Titan-Pinzetten

**ES** Pinzas de titanio


Page 115 to 116



ATRAUMA



TITANIUM



OI 10-2401-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2402-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2403-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2404-00  
Vena  
Closing pressure  
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



OI 10-2405-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2406-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2407-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2408-00  
Vena  
Closing pressure  
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



OI 10-2409-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2410-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2411-00  
Vena  
Closing pressure  
110 +/- 10 gr

OI 10-2412-00  
Vena  
Closing pressure  
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



ATRAUMA



DEBAKEY



OI 10-2413-16  
 16 cm/6 1/4"  
 1.5 mm



OI 10-2413-20  
 20.5 cm/8 1/8"  
 1.5 mm



OI 10-2413-24  
 24.5 cm/9 5/8"  
 1.5 mm



OI 10-2414-16  
 16 cm/6 1/4"  
 1.8 mm



OI 10-2414-20  
 20.5 cm/8 1/8"  
 1.8 mm



OI 10-2414-24  
 24.5 cm/9 5/8"  
 1.8 mm

